










Size & Weight Range Guide (approximately)

Please note the following sling guide is a recommendation only. A full risk assessment must be done prior to any sling being selected. This will ensure safety for the patient and caregiver.

| Size | Low | High | Height |
|------------------|-----|------|---------|
| Small (S) | 75 | 150 | 59"-64" |
| Medium (M) | 125 | 200 | 63"-68" |
| Large (L) | 175 | 300 | 67"-72" |
| Extra Large (XL) | 275 | 500 | 71"-76" |








AdvaCare Systems Stocked Slings

| | | | |
|---|---|-----------------|---|
| <p>Standard Sling Full Electric Lift (Full Body)</p> |  <p>Full Electric Lift</p> | <p>M, L, XL</p> |  |
| <p>Commode Sling Full Electric Lift (Full Body)</p> |  <p>Full Electric Lift</p> | <p>M, L, XL</p> |  |
| <p>Standing Sling</p> |  <p>Sit-to-Stand</p> | <p>M, L, XL</p> |  |
| <p>Transport Sling (Transfer Sling)</p> |  <p>Sit-to-Stand</p> | <p>M, L, XL</p> |  |

Additional slings and sizes can be ordered as a specialty purchase

Sit-to-Stand



Electric Lift

